



Push Down Table Top Scissors

CTT-1

These scissors are easy to use simply requiring a pushing down action on the "T" shaped handle. The spring will automatically reopen the scissors once the pressure is released.

Can be used either by maintaining the paper in a fixed position and "scooting" the scissors along with each cutting action; or by keeping the scissors stable (a Dycem® mat can help with this) and moving the paper towards the scissors.



NOW SUPPLIED IN HANDY ZIP CASE



OUR OT FINDS THESE GOOD FOR...

- ✓ Weak hand and/or arm strength
- ✓ Poor motor control
- ✓ Tremor
- ✓ Joint protection in conditions like Arthritis
- ✓ Cost effective table top scissor that can be stored in regular classroom storage space
- ✓ Children with poor strength able to access/carry the scissor independently

FEATURE	BENEFITS
NEW! Hidden spring	This hidden spring automatically reopens the scissors after every cut
Large "T" shaped handles	Designed for use on table-top to give good stability and leverage
High quality stainless steel blade	Provides longevity and cuts through paper with ease
Blade guard	Fixes over the blades to provide safe storage
NEW! Zip case	Attractive and handy for storing other small items such as pencil grips. Keeps the scissors safe when not in use.

RECOMMENDED FOR



Weak Grip



Poor Hand Control



Tremor